OPERATING INSTRUCTIONS

EN

1. Twenty to thirty minutes before use, remove the product box from the refrigerator, remove the blister pack from the box and allow the syringe to come to room temperature. Do not use if the sterile package appears to be damaged. Use aseptic technique when handling the syringe and administering the injection.

2. **Peel of the blister backing.** (Fig. 1)

3. While holding the blister open side down, bend the blister and allow the syringe to fall gently onto the sterile surface. (Fig. 2)

Alternatively, hold the blister open side up and bend back the blister until the barrel's luer end is exposed. Gripping the luer end of the barrel, remove the syringe from the blister. **Do not remove the syringe from the plunger end.**

4. **Remove the tip cap** (Fig. 3.1, Fig. 3.2) **from the syringe and attach the cannula** (Fig 4.1, Fig. 4.2). **Attention:** Do not apply pressure to the plunger rod while the cannula is being affixed. Verify that the cannula is properly locked to the Luer Lock Adaptor (LLA). Do not overtighten the LLA; this can lead to loosening of the LLA from the barrel.

5. **Apply gentle pressure to the plunger** in order to expel air from the cannula and verify that the syringe is operating properly. (Fig. 5)

6. The syringe is ready for use.

7. Do not reuse or resterilize.

Safe disposal of device : The syringe and cannula should be disposed of as is done for other similar medical products. Used syringes may be disposed of with regular trash or as per glass waste. Cannulas may be disposed of with regular trash or as per sharps waste.

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
6309501203

BioLon®

	Caution, see instructions for use	LOT	Batch code
	See instructions for use		Use by (YYYY-MM-DD : year-month-day)
	Do not reuse	REF	Catalog number
	Protect from light	MD	Medical device
	Do not use if the packaging has been opened or damaged		Keep dry
STERILE EO	Sterilized by ethylene oxide	2°C 8°C	Temperature limit
STERILE A	Sterilized by aseptic techniques		Do not resterilize
	Manufacturer		Double sterile barrier system

MANUFACTURED BY
Nordic Group B.V.
Siriusdreef 41,
2132 WT Hoofddorp,
The Netherlands

DISTRIBUTED BY
NORDIC
PHARMA

6309501203 Issue date: 06/2025 - v1

STERILE EO

CE 2797

BioLon®

BIOLON®	EN										
Sodium hyaluronate, 10 mg/ml	PRODUCT INFORMATION										
FORMULA Each ml of BIOLON® contains: <table> <tr> <td>Sodium hyaluronate</td><td>10 mg</td></tr> <tr> <td>Sodium chloride</td><td>8,5 mg</td></tr> <tr> <td>Disodium hydrogen phosphate dodecahydrate</td><td>0,56 mg</td></tr> <tr> <td>Sodium dihydrogen phosphate dihydrate</td><td>0,05 mg</td></tr> <tr> <td>Water for injections</td><td>q.s.</td></tr> </table>	Sodium hyaluronate	10 mg	Sodium chloride	8,5 mg	Disodium hydrogen phosphate dodecahydrate	0,56 mg	Sodium dihydrogen phosphate dihydrate	0,05 mg	Water for injections	q.s.	
Sodium hyaluronate	10 mg										
Sodium chloride	8,5 mg										
Disodium hydrogen phosphate dodecahydrate	0,56 mg										
Sodium dihydrogen phosphate dihydrate	0,05 mg										
Water for injections	q.s.										

DESCRIPTION

Sodium hyaluronate is a physiological material found in both animal and human tissues. It is a polysaccharide consisting of a repeating disaccharide of N-acetylglucosamine and sodium glucuronate, linked by alternating β-1,3 and β-1,4 glycosidic bonds. Chemically identical in all species, hyaluronate can be found in the vitreous and aqueous humor of the eye, in synovial fluid, in skin and in the umbilical cord. BIOLON® is a sterile 1% (10 mg/ml) solution of highly purified, high molecular weight viscoelastic sodium hyaluronate in phosphate-buffered saline (pH 6,8-7,6, Osmolality 258-381 mOsm/kg). When introduced into the anterior segment of the eye during surgical procedures, BIOLON® maintains a deep anterior chamber throughout surgery, facilitates surgery and reduces trauma to the corneal endothelium and surrounding tissues. Its viscoelasticity also helps to repel the vitreous face and discourage formation of a postoperative flat chamber. In posterior segment surgery, BIOLON® serves as a surgical aid in separating, moving, and holding tissues. BIOLON®creates a clear field of vision for inspection of the retina and for photocoagulation during and after surgery. It does not interfere with epithelialisation and normal wound healing. Any traces of BIOLON® left in the anterior segment of the eye after surgery dissipate within a week.

INDICATION

Eye surgery, including intraocular lens insertion, intracapsular and extracapsular lens extraction, glaucoma surgery, corneal graft surgery for accidental trauma, retinal detachment and vitreal replacement procedures.

CONTRAINDICATIONS

When used as recommended, there are no known contraindications.

PRECAUTIONS

Precautions are limited to those normally associated with the surgical procedures being performed.Instilling excessive amounts of BIOLON® into the anterior or posterior segments of the eye may increase intraocular pressure. Increased postoperative intraocular pressure may also be caused by a pre-existing glaucoma condition or by compromised outflow and by operative procedures and sequelae thereof, including enzymatic colonysis, absence of an iridectomy, trauma to filtration structures, and by blood and lenticular remnants in the anterior chamber. Because these factors vary from case to case and are difficult to predict, the following precautions are recommended:

Do not overfill the eye with BIOLON®. At the end of surgery, all remaining BIOLON® should be removed by irrigation or aspiration. Carefully monitor the intraocular pressure, especially during the immediate postoperative period. If a significant increase in pressure is observed, treat appropriately. Avoid the trapping of air bubbles behind BIOLON®. Since the presence of minute quantities of impurities (proteins, etc.) cannot be totally excluded, the physician should be aware of immunological, allergic and other potential risks associated with the injection of biological substances. On rare occasions, viscoelastic products containing sodium hyaluronate have been known to become slightly opaque or to form a slight precipitate upon instillation into the eye. The clinical significance, if any, of this phenomenon is not known. The physician should, however, be aware of this possibility, and any visible particulate or precipitated material should be removed by irrigation and/or aspiration.

Use of the cannula: the product is for single use. Reuse of single use device can result in injury or impairment requiring professional medical intervention as the safety and effectiveness cannot be guaranteed.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DOSAGE

Dosage varies with type of surgery. Usually a dose of 0.2 to 0.6 ml is injected into the anterior segment of the eye. Greater amounts are used in the posterior segment.

EXAMPLES OF SURGICAL APPLICATIONS

Cataract Surgery: Before lens extraction, BIOLON® is introduced into the anterior chamber to protect the corneal endothelium and to maintain a deep anterior chamber. Additional amounts may also be introduced before implantation of an intraocular lens, and to coat the artificial lens and the surgical instruments.

Glaucoma Surgery: BIOLON® is instilled before the trabeculectomy in order to reconstitute the anterior chamber. At the close of surgery, further instillation of BIOLON® may be required to achieve good sub- conjunctival filtration and to prevent tissue adhesences.

Corneal Transplant Surgery: In trepanation of the cornea, the anterior chamber is filled with BIOLON®. The donor graft should be placed on the surface of the solution and sutured into position. Additional amounts BIOLON® may be injected to maintain a deep anterior chamber. BIOLON® may also be used in the anterior chamber of the donor eye to protect the corneal endothelial cells of the graft prior to trepanation, and to protect the exposed endothelial layer of the donor button during preparation of the recipient eye.

Retinal and Vitreous Surgery: BIOLON® should be introduced slowly into the vitreous cavity. By carefully directing the injection, BIOLON® can be used to separate membranes from the retina for safe excision and release of traction. BIOLON® also helps in moving tissues into desired positions. For example, it is helpful in gently pushing back a detached retina, in unrolling a retinal flap, or in holding the retina against the sclera during reattachment.

INCOMPATIBILITIES

Mixing of quaternary ammonium salts such as benzalkonium chloride with sodium hyaluronate solutions results in formation of a precipitate. BIOLON® should not be administered through a cannula previously used with medical solutions containing benzalkonium chloride.

ADVERSE REACTIONS

BIOLON® is well-tolerated in the human eye. Transient rises of post-operative intraocular pressure have been reported in some cases. A causal relationship has not been established between BIOLON® use and postoperative inflammatory reactions (iritis), corneal edema and corneal decompensation.

INTERACTIONS

None currently known.

PRESENTATIONS

BIOLON® is supplied in sterile disposable syringes containing 0,5 ml or 1.0 ml solution.

SHELF LIFE

3 years. See product label for expiration date.

STORAGE INSTRUCTIONS

Store at 2°- 8°C. May be kept at 25°C for up to six months. Do not freeze. Protect from light.

To identify, access and view the electronic Instructions For Use (IFU), go to the following URL: https://www.nordicpharma.com/our-products/ and access the BIOLON section

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